










## SYMBOLS GLOSSARY FOR ASTROLABE

*This document contains symbols that may appear on the product packaging or labeling. It provides information on the symbol including the standards title, reference, and description, where available. This document does not contain company created symbols, pictograms, or graphical images.*

*Please read the full Instructions For Use for detailed information on proper use, indications, contraindications, warnings, precautions and adverse events. This document does not replace the Instructions For Use.*

**ISO 15223-1 Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements**

Symbol	Description	Chapter	Description
	MRI Conditional	5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Batch Code	5.1.5	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Caution	5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	Consult Instructions of Use	5.4.3	Indicates the need for the user to consult the instructions for use.
	Manufacturer	5.1.1	Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. NOTE 5 The date of manufacture, as well as the name and address of the manufacturer, can be combined in one symbol.
	Do not re-use	5.4.2	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	Non-Sterile	5.2.7	Indicates a medical device that has not been subjected to a sterilization process.
	Do not use if package is damaged	5.2.8	Indicates a medical device that should not be used if the package has been damaged or opened.
	Date of Manufacturer	5.1.3	Indicates the date when the medical device was manufactured.