

7.0mm

CANNULATED SCREW | TORNILLO CANULADO

Box and Trays
Caja y Estuches

Ref.	Box and Trays Caja y Estuches
09.50.03.00130	Box w/ Brackets, 250 x 330 x 160mm, Cannulated Screw 7.0mm <i>Caja con Soportes, Tornillo Canulado 7.0mm</i>
09.50.04.00210	Implant Case, Washers 7.0mm <i>Estuche para Implantes, Arandelas 7.0mm</i>
09.50.06.00130	Instrument Tray, Set-1, Cannulated Screw 7.0mm <i>Bandeja Instrum, Set-1, Tornillo Canulado 7.0mm</i>
09.50.06.00131	Instrument Tray, Set-2, Cannulated Screw 7.0mm <i>Bandeja Instrum, Set-2, Tornillo Canulado 7.0m</i>
09.50.10.00230	Dispenser for K-Wire 230mm <i>Dispensador para Agujas de Kirschner 230mm</i>

Cannulated Screw, Fully Threaded
Tornillo Canulado, Rosca Total

Ref.	Length Longitud
15.30.01.70020	20mm
15.30.01.70025	25mm
15.30.01.70030	30mm
15.30.01.70035	35mm
15.30.01.70040	40mm
15.30.01.70045	45mm
15.30.01.70050	50mm
15.30.01.70055	55mm
15.30.01.70060	60mm
15.30.01.70065	65mm
15.30.01.70070	70mm
15.30.01.70075	75mm
15.30.01.70080	80mm
15.30.01.70085	85mm
15.30.01.70090	90mm
15.30.01.70095	95mm
15.30.01.70100	100mm
15.30.01.70105	105mm
15.30.01.70110	110mm
15.30.01.70115	115mm
15.30.01.70120	120mm
15.30.01.70125	125mm
15.30.01.70130	130mm



7.0mm

CANNULATED SCREW | TORNILLO CANULADO

Cannulated Screw, 16mm Threaded
Tornillo Canulado, Rosca 16mm

Ref.	Length <i>Longitud</i>
15.44.01.70030	30mm
15.44.01.70035	35mm
15.44.01.70040	40mm
15.44.01.70045	45mm
15.44.01.70050	50mm
15.44.01.70055	55mm
15.44.01.70060	60mm
15.44.01.70065	65mm
15.44.01.70070	70mm
15.44.01.70075	75mm
15.44.01.70080	80mm
15.44.01.70085	85mm
15.44.01.70090	90mm
15.44.01.70095	95mm
15.44.01.70100	100mm
15.44.01.70105	105mm
15.44.01.70110	110mm
15.44.01.70115	115mm
15.44.01.70120	120mm
15.44.01.70125	125mm
15.44.01.70130	130mm
15.44.01.70135	135mm
15.44.01.70140	140mm
15.44.01.70145	145mm
15.44.01.70150	150mm



7.0mm CANNULATED SCREW | *TORNILLO CANULADO*

Cannulated Screw, 32mm Threaded *Tornillo Canulado, Rosca 32mm*

Ref.	Length <i>Longitud</i>
15.46.01.70045	45mm
15.46.01.70050	50mm
15.46.01.70055	55mm
15.46.01.70060	60mm
15.46.01.70065	65mm
15.46.01.70070	70mm
15.46.01.70075	75mm
15.46.01.70080	80mm
15.46.01.70085	85mm
15.46.01.70090	90mm
15.46.01.70095	95mm
15.46.01.70100	100mm
15.46.01.70105	105mm
15.46.01.70110	110mm
15.46.01.70115	115mm
15.46.01.70120	120mm
15.46.01.70125	125mm
15.46.01.70130	130mm
15.46.01.70135	135mm
15.46.01.70140	140mm
15.46.01.70145	145mm
15.46.01.70150	150mm



7.0mm

CANNULATED SCREW | TORNILLO CANULADO



INSTRUMENTS INSTRUMENTOS		
09.01.01.20100	Drill Bit, Ø2.0 x 100mm, twist length 30mm, AO Coupling <i>Broca Ø2.0 x 100mm, H-30, AO</i>	2
09.01.05.45240	Drill Bit, Cannulated, Øext. 4.5 x Øint. 2.1 x 240mm, twist length 40mm, SZ Coupling <i>Broca Canulada Ø4.5 x 240mm, H-40, SZ</i>	2
09.02.03.70000	Countersink, Cannulated, Ø7.0mm, SZ Coupling <i>Avellanador Canulado Ø7.0mm, SZ</i>	1
09.03.06.70000	Tap, Cannulated, Ø7.0mm, SZ Coupling <i>Macho Canulado Ø7.0mm, SZ</i>	1
09.04.08.00050	T-Handle, Cannulated, SZ Coupling, Black <i>Mango en T Canulado, SZ, Neg</i>	1
09.05.01.20010	Drill Guide, Ø2.0mm, Parallel - Multiple <i>Guía de Perforación Ø2.0mm, Paralelo - Múltiple</i>	1
09.05.01.45000	Drill Guide, Ø4.5mm, Parallel <i>Guía de Perforación Ø4.5mm, Paralelo</i>	1
09.05.02.45320	Drill Guide, Ø4.5/Ø3.2mm, Double <i>Guía de Perforación Ø4.5/Ø3.2mm, Doble</i>	1
09.06.00.45099	Drill Sleeve, Ø4.5 x 99mm <i>Manga de Perforación Ø4.5 x 99mm</i>	1
09.06.02.20108	Wire Sleeve, Ø2.0 x 108mm <i>Manga Ø2.0 x 108mm, para Aguja</i>	2
09.06.04.45000	Holding Sleeve, Ø4.5mm and Larger <i>Manga de Sujeción Ø4.5mm y Superiores</i>	1
09.06.06.11080	Protection Sleeve, Ø11.0/8.0mm <i>Manga de Protección Ø11.0/8.0mm</i>	1
09.07.12.25200	Shaft Screwdriver, Cannulated, Torx-25, 200mm, SZ Coupling <i>Destornillador Canulado, TX-25, 200mm, SZ</i>	2
09.08.02.20160	Gauge 160mm, Direct Measuring, for Wire 230mm, Ø up to 2.0mm <i>Medidor 160mm, Aguja Ø2.0 x 230mm</i>	1
09.10.01.00010	Screw Forceps, Large <i>Pinza para Tornillos, Grande</i>	1
09.25.01.01200	Cleaning Wire, Ø2.0mm, Stylet <i>Alambre de Limpieza Ø2.0mm, Estilete</i>	1
09.25.01.02210	Cleaning Wire, Ø2.1mm, Brush <i>Alambre de Limpieza Ø2.1mm, Cepillo</i>	1
09.26.00.20118	Trocar, Ø2.0 x 118mm <i>Trocar Ø2.0 x 118mm</i>	2

Ø7.0MM CANNULATED SCREWS, FULLY THREADED Ø7.0MM TORNILLO CANULADO, ROSCA TOTAL		
15.30.01.70020	Cannulated Screw, Fully Threaded Ti, Ø7.0 x 20mm Torx <i>Tornillo Canul, Ti, Rosca Total, Ø7,0 x 20mm TX</i>	1
15.30.01.70025	Cannulated Screw, Fully Threaded Ti, Ø7.0 x 25mm Torx <i>Tornillo Canul, Ti, Rosca Total, Ø7,0 x 25mm TX</i>	1
15.30.01.70030	Cannulated Screw, Fully Threaded Ti, Ø7.0 x 30mm Torx <i>Tornillo Canul, Ti, Rosca Total, Ø7,0 x 30mm TX</i>	1
15.30.01.70035	Cannulated Screw, Fully Threaded Ti, Ø7.0 x 35mm Torx <i>Tornillo Canul, Ti, Rosca Total, Ø7,0 x 35mm TX</i>	1
15.30.01.70040	Cannulated Screw, Fully Threaded Ti, Ø7.0 x 40mm Torx <i>Tornillo Canul, Ti, Rosca Total, Ø7,0 x 40mm TX</i>	1

